Novel Coronavirus (2019-nCoV) Antigen Test Kit (Colloidal Gold Immunochromatography)

- Nasal
English

【PRODUCT NAME】

Common name: Novel Coronavirus (2019-nCoV) Antigen Test Kit (Colloidal Gold Immunochromatography)

【SPECIFICATIONS】

1 Test/Kit, 5 Tests/Kit, 20 Tests/Kit

【INTENDED USE】

The COVID-19 Antigen Rapid Test Kit is intended for in vitro qualitative detection of SARS-CoV-2 nucleocapsid antigens in nasal swab from individuals who are suspected of COVID-19.

The novel coronaviruses (SARS-CoV-2) belong to the β genus. COVID-19 is an acute respiratory infectious disease. Based on the current epidemiological investigation, the incubation period is 1 to 14 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions.

Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19, a confirmed with a molecular assay, if necessary, for patient management.

【PRINCIPLE】

The COVID-19 Antigen Rapid Test Kit is a lateral flow immunoassay based on the principle of the double-antibody sandwich technique. The test uses SARS-CoV-2 antibody (test line T) and goat anti-mouse IgG (control line C) immobilized on a nitrocellulose strip. The burgundy-colored conjugate pad contains colloidal gold conjugated to anti-SARS-CoV-2 antibody conjugated with colloidal gold (SARS-CoV-2 conjugates) and mouse IgG-gold conjugates.

During the test, SARS-CoV-2 antigen in the specimen interacts with SARS-CoV-2 antibody conjugated with color microparticles making antigen-antibody labeled complex. This complex migrates on the membrane via capillary action until the test line, where it will be captured by the pre-coated SARS-CoV-2 nucleocapsid protein monoclonal antibody. A colored test line (T) will be visible in the result window if COVID-19 antigens are present in the specimen. Absence of the T line suggests a negative result. The control line (C) is used for procedural control and should always appear if the test procedure is performed properly.

【MATERIALS PROVIDED】

NO.	Composition	For 1 Test/Kit	For 5 Tests/Kit	For 20 Tests/Kit
1	2019-nCoV Antigen Detection Card (sealed foil pouch)	1	5	20
2	Sample extraction buffer	1	5	20
3	Sterile swab	1	5	20
4	Tube stand	-	1	1

Note: 1. Reagents from different batches cannot be mixed use.

【STORAGE AND STABILITY】

- The storage condition of this kit is 2-30°C in a dry and dark place. The shelf life is 24 months.
- Once open the pouch, the test should be used within one hour.
- The LOT at the expiration date were printed on the labeling.

【SPECIMEN】

Nasal swab specimen collection

Take the swab out of the package and carefully insert it into the nasal cavity 1.5cm until you feel a slight resistance. If you feel strong resistance or pain, do not insert the swab deeper. Under moderate pressure, move the swab in a circular motion along the inner wall of the nose for at least 15 seconds (approximately 4-6 times) to collect as many cells and mucus as possible. Repeat the sample with the same swab in the other nostril.

Sample transport and storage

The virus sampling solution or the sample extraction solution provided by this kit should be used for processing as soon as possible after specimen collection. If it cannot be processed immediately, the specimen should be stored in a dry, sterilized and strictly sealed plastic tube immediately. It can be stored for 8 hours at 2°C-8°C or -20°C can be stored for a long time.

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1. Preparation

Allow the test devices, reagents to equilibrate to room temperature (15-30°C) prior to testing.


2. Extraction

- Tear off the sealing film of the extraction tube.
- Insert the sampled swab into the solution in the sample extraction tube, rotate it close to the inner wall of the test tube about 10 times, Leave the swab in extraction buffer for 1 minute.


- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- Install the cap onto the sample extraction tube tightly.

3. Test

- Open the pouch and remove the card. Place the card on a flat and level surface.
- Add 4 drops (about 100 μ l) of the processed sample extract to the sample hole of the test card.
- Read the test results after 15-20 minutes. Results after 20 minutes may not be accurate.

【INTERPRETATION OF RESULTS】

Positive

If two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive. No matter how faint the colored band is in the Test Zone (T), the result should be considered as positive.

Negative

If one colored band appears in the Control Zone (C) and no colored band appears in the Test Zone (T) within 15-20 minutes, the test result is negative.

Invalid

If no color line appears in the control area (C) within 15-20 minutes, the test is invalid.

【LIMITATIONS】

- This product is only used to test the secretions of individual nasal swabs.
- The accuracy depends on the sample collection process. Improper sample collection, storage, repeatedly frozen and thawed will affect the result.
- The result of this product is for clinical reference only, and should not be used as the only basis for clinical diagnosis and treatment. The information of patient's symptoms, signs, history, other tests (especially pathogenic tests), treatment response and epidemiology are necessary to comprehensive consideration.
- The presence of drugs in the samples such as high concentrations of nasal sprays will interfere with the results. If the result is suspicious, please retest.
- Weak positive samples that have been heat-inactivated and treated with guanidine salt may cause false negative results.
- The possibility of false negative results:
 - Improper operation procedures (sample collection, transfer, storage, processing, etc.) or the virus in the sample is below the detection limits.
 - Possible subtypes (mutant strains).
 - The antigen that are detected are not present during the stage of disease in which a sample is collected.

【PERFORMANCE CHARACTERISTICS】

1. When testing with enterprise references, meet the following standards:

- Negative references compliance rate: Use the enterprise negative references for testing, and the negative references should be detected at least 20/20 (-/-).
- Positive references compliance rate: Use the enterprise positive references for testing, and the positive references should be detected at least 5/5 (+/+).
- Sensitivity references: Each detection item is tested with the corresponding minimum detection limit quality control serum. L1 and L2 results should be positive, and L3 should be negative.
- Repeatability: Use enterprise precision references for testing, and the test results of repeatable references should be consistent.

2. Limit of detection

The study used cultured SARS-CoV-2 virus, which is β -propiolactone and heat inactivated and spiked into sample extraction buffer. The limit of detection is 9.65 TCID₅₀/ml.

3. High-dose hook effect

The Novel Coronavirus (2019-nCoV) Antigen Test Kit was tested up to 7 \times 10⁷ TCID₅₀/ml of inactivated SARS-CoV-2 and no high-dose hook effect was observed.

4. Cross reactivity

Cross reactivity was evaluated by testing 56 commensal and pathogenic microorganisms that may be present in the throat. Each of the organisms and virus were tested in the absence or presence of heat inactivated SARS-CoV-2 virus (26S TCID₅₀/ml), no cross-reactivity was seen.

No.	Crossing reacting substance	Strain	Concentration of cross-reacting substance	Result	
				Absence of SARS-CoV-2	Presence of SARS-CoV-2
1	Human Coronavirus	HKU1	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
2		229E	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
3		OC43	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
4		NL63	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
5	SARS	SARS	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
6		MERS	2 \times 10 ⁷ TCID ₅₀ /ml	-	+

7	Adenovirus	Type 1	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
8		Type 2	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
9		Type 3	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
10		Type 4	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
11		Type 5	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
12		Type 7	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
13		Type 55	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
14	Human Metapneumovirus (hMPV)	hMPV 3 Type B1/Peru2-2002	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
15	Parainfluenza virus	hMPV 16 Type A1/A10-2003	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
16		Type 1	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
17		Type 2	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
18	Influenza A	Type 3	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
19		Type 4A	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
20		H1N1	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
21		H3N2	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
22	Influenza B	H5N1	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
23		H7N9	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
24		Yamagata	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
25	Enterovirus	Victoria	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
26		Type 68	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
27	Respiratory syncytial virus	09/2014 isolate 4	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
28		Type A	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
29	Rhinovirus	Type B	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
30		A16	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
31	Chlamydia pneumoniae	Type B42	2 \times 10 ⁷ TCID ₅₀ /ml	-	+
32		TWAR strain TW-183	5 \times 10 ⁷ CFU/ml	-	+
33	Haemophilus influenzae	NCTC 4560	5 \times 10 ⁷ CFU/ml	-	+
34		Bloomington-2	5 \times 10 ⁷ CFU/ml	-	+
35	Legionella pneumophila	Los Angeles-1	5 \times 10 ⁷ CFU/ml	-	+
36		82A3105	5 \times 10 ⁷ CFU/ml	-	+
37	Mycobacterium tuberculosis	K	5 \times 10 ⁷ CFU/ml	-	+
38		Erdman	5 \times 10 ⁷ CFU/ml	-	+
39		H37Rv	5 \times 10 ⁷ CFU/ml	-	+
40		CDC1551	5 \times 10 ⁷ CFU/ml	-	+
41	Streptococcus pneumoniae	H37Rv	5 \times 10 ⁷ CFU/ml	-	+
42		4752-	5 \times 10 ⁷ CFU/ml	-	+
43	Streptococcus pneumoniae	98[Maryland(01)68-17]	5 \times 10 ⁷ CFU/ml	-	+
44		178[Poland 23F-16]	5 \times 10 ⁷ CFU/ml	-	+
45		262[CIP 104340]	5 \times 10 ⁷ CFU/ml	-	+
46	Streptococcus pyogenes	Slovakia 14-1(20)9055]	5 \times 10 ⁷ CFU/ml	-	+
47		Typing strain T1(NCIB13841, SF 130)	5 \times 10 ⁷ CFU/ml	-	+
48	Bordetella pertussis	NCCP 13671	5 \times 10 ⁷ CFU/ml	-	+
49		Mutant 22	5 \times 10 ⁷ CFU/ml	-	+
50	Mycoplasma pneumoniae	FH strain of Eaton Agent (NCTC 10119)	5 \times 10 ⁷ CFU/ml	-	+
51		M129-87	5 \times 10 ⁷ CFU/ml	-	+
52	Pneumocystis jirovecii (PJP)	N/A	N/A	-	+
53	Pooled human nasal wash	N/A	N/A	-	+
54	Candida albicans	3147	5 \times 10 ⁷ CFU/ml	-	+
55	Pseudomonas aeruginosa	R. Hugh 813	5 \times 10 ⁷ CFU/ml	-	+
56	Staphylococcus epidermidis	FDA strain PC1 1200	5 \times 10 ⁷ CFU/ml	-	+
57	Streptococcus salivarius	5218 [FO 13956]	5 \times 10 ⁷ CFU/ml	-	+

5. Interference

There was no interference for potential interfering substances listed below.

1) Exogenous factor

No	Exogenous factor	Interfering substance	Test conc
1	Nasal sprays or drops	Phenylephrine	128 μ g/ml
2		Oxymetazoline	128 μ g/ml
3		Saline Nasal Spray 10%	10%(V/V)
4		Dexamethasone	2 μ g/ml
5	Nasal corticosteroids	Flunisolide	0.2 μ g/ml
6		Triamcinolone acetonide	0.2 μ g/ml
7		Mometasone	0.3 μ g/ml
8	Throat lozenges	Strepsils (flurbiprofen 8.75mg)	5% (w/v, 50mg/ml)
9		Throat candy	5% (w/v, 50mg/ml)
10	Oral anesthetic	Amesol (Benzocaine 20%)	5%(V/V)
11		α -interferon-2b	0.01 μ g/ml
12		Zanamivir (influenza)	2 μ g/ml
13		Ribavirin (HCV)	0.3 μ g/ml
14		Oseltamivir (influenza)	2 μ g/ml
15		Peramivir (influenza)	60 μ g/ml
16		Lopinavir (HIV)	80 μ g/ml
17		Ritonavir	20 μ g/ml
18		Arbidol (influenza)	40 μ g/ml

19	Antibiotic	Levofloxacin Tablets	40 μ g/ml
20		Azithromycin	200 μ g/ml
21		Ceftriaxone	800 μ g/ml
22		Meropenem	100 μ g/ml
23		Tobramycin	128 μ g/ml
24	Other	Mucin: bovine submaxillary gland, type	100 μ g/ml
25		Biotin	100 μ g/ml
26	antihistamine drug	Diphenhydramine Hydrochloride	2 μ g/ml
27		Loratadine	40 μ g/ml

2) Endogenous factor

No	Exogenous factor	Interfering substance	Test conc
1	Autimmune disease	Human anti-mouse antibody, HAMA	800 μ g/ml
2	Serum protein	Whole Blood(human), EDTA, anticoagulated	10%(w/w)

6. Clinical performance

Nasal swab specimen

To estimate the clinical performance between the Novel Coronavirus (2019-nCoV) Antigen Test Kit and the RT-PCR comparator, 210 samples were collected. The result is below:

COVID-19 Antigen	RT-PCR		Total
	Positive	Negative	
Positive	102	3	103
Negative	6	101	107
Total	108	102	210

Sensitivity (PPA)=94.44% (102/108), (95%CI:90.12%-98.76%)

Specificity (NPA)=99.02% (101/102), (95%CI: 97.11%-100%)

Total consistency=96.67% (203/210), (95%CI:94.24%-99.09%)

【WARNINGS AND PRECAUTIONS】

- For in vitro diagnostic use only.
- Please read the instruction manual carefully and strictly follow the methods and steps to perform the test.
- This product is a qualitative detection; results cannot be used as a quantitative basis.
- Do not use beyond the expiration date.
- Test is for single use only. Do not re-use under any circumstances.
- Due to the different sample tier, the red test line will show different shades of color, which are all represented positive results. The test line color shades are not applicable for determining the level of antigen tier in the sample.
- Cold temperature preserved sample should be equilibrated to room temperature and mix thoroughly before testing.
- Pay attention to safety measures during operation, such as protective clothing and gloves.

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	Keep dry		Keep away from sunlight		In vitro diagnostic use
	Manufacturer		Authorized representative		Date of manufacture
	Batch number		Store between 2-30°C		Expiry Date
	Do not use if package is damaged		See instructions for use		Catalogue number
	Do not re-use		Contains sufficient for <n> tests		CE mark

【BASIC INFORMATION】

Gallergo Co., Ltd.
Address: North side of the third floor of No. 6 factory building, No. 1 Xining Road, Wuqing Development Zone, Tianjin
Contact: 86-22-82939117

Omunda Medical Technology Service GmbH
Address: Treskowallee 108, 10318 Berlin, Germany
Contact: 0049-30-81865123

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